DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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emergency.

- ◆ Never try to repair the Thermal Cycler by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the Thermal Cycler (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the equipment and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.
- ◆ Clean the block's wells whenever anything is spilled into them or dust gathers.
- ◆ Clean the wells with a swab moistened with water and 95% ethanol. If using sodium hypochlorite, swab the wells with water afterward to remove all traces of it. Do not clean the block with caustic or strongly alkaline solutions (e.g., strong soaps, ammonia, sodium hypochlorite) at a higher concentration. These can damage the block's protective anodized coating.

Note: Be sure to shut-off the Thermal Cycler when cleaning the surface and the block.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INSTRUCTIONS ON ENVIRONMENT PROTECTION

Do not dispose of this product in the usual household garbage at the end of its life cycle; hand it over at a collection point for the recycling of electrical and electronic appliances. The symbol on the product, the instructions for use or the packing will inform about the methods for disposal.

The materials are recyclable as mentioned in its marking. By recycling, material recycling or other forms or re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location

For example, if you hope the lid turn off when the sample is hold at 4 °C, press [4] and [ENTER]. The value will not lost unless it is modified again.

If you are not using the Thermal Cycler for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

Security

- ◆ The equipment must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ You should put the Thermal Cycler in a horizontal plane stable table, having a safety area of at least 30 cm per side.
- ◆ Do not place the Thermal Cycler near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ During its functioning dangerous materials such as flammable or pathological substances must be out of the safety area.
- ◆ When you are not using the Thermal Cycler for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the Thermal Cycler by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the Thermal Cycler (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The equipment is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of

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


Press [▲] or [▼] to page up and page down fast to the pages you want to modify, press [CANCEL] to clear the old value before you input a new one, and must press [ENTER] to accept a new data. Lastly the LCD window will display:

If “Yes” is selected, the old program will be replaced by the new one after several minutes.


If “No” is selected, modified program will be lost. The old program will be reserved.

File


To delete a stored program, select **File** from the main menu and press [ENTER], then the LCD window will show:

Press [ENTER] and all programs stored will be listed.


Select the program you want to deleted, press [ENTER]. If select “Yes”, the program will be deleted. If select “No”, the program will be reserved.

Lid

To select a minimum block temperature below which the heated **lid** will automatically turn off, select Lid from main menu, then press [ENTER], the following screen will be displayed.

Note: When Editing a new program, you might press [STOP] accidentally, the screen will ask you if save it, just press [▲] or [▼] to return to the former page.

At this time, a new program is finished. From main menu you could select List or Edit to review or modify it lately.

List

To view a program in the LCD window, select List from the main menu, then press [ENTER], the screen will show all the program stored:

LIST:	
▲1	a2
a3	a4
a5	a6

Locate the program to be viewed, then press [ENTER], the program will immediately be listed in the LCD window. Listed programs cannot be modified.

The last page will display:

LIST:	a1
Estimated run time:	
99h99m99s	
Save? yes	NO

Only [No] could be selected at this time, press [ENTER] to return to the main menu.

To modify programs, select *Edit* from the main menu. To delete programs, select *File* from the main menu

Edit

To initiate editing, select Edit from the main menu, then press [ENTER]. The screen will display:

EDIT:	
▲1	a2
a3	a4
a5	a6

Select the program to be modified, then press [ENTER]. The first editing screen display control method first.

EDIT:	a1
Control Method:	
BLOCK	sim-

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
Hold at 72.0C
Hold for 10m 0s

Press [7] [2] [0] and [ENTER] and then press [1] [0] and [ENTER]

Press [ENTER] and the screen will display:

Segment # 4
hold cycle END

Hold at 0.0 C
Hold for 0m 0s

If you need 4 °C hold for forever, move the cursor to **Hold** and press [ENTER]

Hold at 4.0 C
Hold for 99m 0s

Press [4] [0] and [ENTER] and press [9] [9]

Note: "Hold for 99 minute" means "Hold forever"

Note:: END must be selected as the final segment of a program

Press enter

Segment # 5
hold cycle END

If you want to finish the program, move the cursor to **End**, press [ENTER] and the screen will say:

Enter: a1
Estimated Run Time:
99h 99m59s
Save: YES no

If Yes is selected, the program will be saved. If No is selected, the program you just enter will lost. The screen will display main menu after several seconds.

MAIN: <4196>
RUN enter
list edit
file lid

Note: The temperature on display represents the middle of the block, so 15 sec added to the normal holding time is needed to generate gradient on block.

If 20 °C gradient temp. is needed, just enter [1] [0] and [ENTER]. Then the cursor will move to the Calculator automatically, press [ENTER] to preview the temp. distribution on block.

Gradient Preview:
Col1 40.0C left →

You could preview the temperature of 12 columns from left to right by pressing [▶] cursor key.

Press [▼] [▲] or [STOP] to return to the below page.

Step # 1
Option? no YES

At this time, you could select No to enter next step.

Step # 2
0.0C Hold 0m0s

As per the above operation, you could manage the next editing. When all steps are set, the screen will display:

Total Cycle = 0

Total Cycle = 35

If 35 cycles is needed, press [3] [5] and

Segment # 3
hold cycle END

Press [ENTER]

Note: You should press [▲] or [▼] to Page Down or Page Up fast

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```
MAIN: <4196>
RUN   enter
lis   edit
files lid
```

Note: Temperature Increment/Decrement could be used in Touchdown PCR experiment. The available range is -10.0 °C - 10.0 °C

If Temp. increment 0.5 °C is needed, press [5]. If Temp. decrement 1 °C is needed, press [-] [1] [0]. If the display is acceptable, press [ENTER].

If Time Increment/Decrement is needed, move the cursor to *Ext*, press [ENTER]

```
Step # 1
Option? no YES
```

If select Time increment 1sec, press [1] and [ENTER]. If Time decrement 1sec is needed, press [-][1] and [ENTER]. If the display is correct, press [ENTER] to return to the option page. Press [▼] [▲] or [STOP], the screen will display like below:

```
Step # 1
Option? no YES
```

If you want the ramping rate at 1 °C/sec to the finish temp. 94 °C, move the cursor to *Ramp*, press [ENTER]

```
Step # 1
Ramp rate = 0.0C/ s
```

Note: Ramp Rate Range is 0.3 °C - 3 °C/ sec. If the 0.0 °C/s -0.2 °C/s is entered, the instrument will run at the max rate. Attention please, the setting point should be the finish temperature instead of start one

Press [1] [0] and [ENTER] and return to the option page.

If you need auto pause function, move the cursor to *Pause* and press [ENTER]

```
Step # 1
Pause? OFF on
```

Note: Auto Pause could be used in hot-start technology.

Select On or Off to confirm and press [▼] or [▲] key to return to the former page.

If select Gradient function (only model 4196), move the cursor to *Gradient*, press [ENTER]

```
Enter: a1
Step#1 Calculator
Gradient = +/- 0 C
```

If you want to select "Cycle", move the cursor to Cycle.

Enter: a1
Segment # 2
hold CYCLE end

Press [ENTER] 0 step PCR

If select 3-step cycle, input [3] 3 step PCR

Press [ENTER] Step # 1
0.0C Hold 0m 0s

If select 94°C hold for 30sec, press [9] [4] [0] and [ENTER], and press [3] [0] and [ENTER], the screen will display:

Step # 1
94.0C Hold 0m 30s

Press [ENTER], and the screen will show:

Step # 1
Option? NO yes

Optional functions:

1. Temp. Increment/ Decrement
2. Time Increment/Decrement
3. Ramp Rate adjustable
4. Auto Pause Settable
5. Gradient Function (Model 4196)

If Option is not needed, move the cursor to No, press [ENTER] to enter the next step.

If Option is needed, move the cursor Yes, press [ENTER] to show an option page

Step # 1
INC ext ramp pause

Enter: a1
Step # 1
Delta temp. = 0.0C

If Temp. Increment/Decrement is selected, move the cursor to Inc, and then [ENTER]

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```
Press<STOP>to abort,
<ENTER> to run again
```

Note: The Screen will show Blk= 94C if block control method is selected, and the screen will show Smp= 94C if sim-tube control method is selected.

When a temperature control method has been chosen, the Enter menu will be displayed:

```
Enter:      a1
Segment # 1
Hold Cycle END
```

If select hold, move the cursor to Hold, then press [Enter]

```
Enter:      a1
Hold at 0.0 C
Hold for 0m 0s
```

If select 94 °C hold for 5 mins, press [9] [4] [0] and press [ENTER] to confirm. At this time, the cursor will move before "m" automatically, press [5] and [ENTER], then the cursor will move before "s"

```
Enter:      a1
Hold at 94.0C
Hold for 5m 0s
```

If acceptable, press [ENTER] and the screen will display:

```
Enter:      a1
Segment # 2
hold cycle END
```

Note: Be sure to press [CANCEL] to clear the old value before you input a new one, and must press [ENTER] to accept a new data.

Note: Temperature available Range: 0 °C - 99.9 °C. Minute available Range: 0min - 99min. Second available Range: 0sec - 59sec.

Note: please shut off the instrument when the screen shows the main menu. Turning off the instrument does not stop a running protocol. The instrument will assume the protocol was stopped by a power outage and will resume running the protocol when the instrument is turned back on.

Warning: when you run the block at 4 °C forever, if the ambient humidity is high, the block's surface might absorb moisture in the air and begin to dew. Under such condition, we suggest change 4 °C to 10 °C or above to avoid condensation water leaking into the cyclor.

Enter

To design a new program, select Enter from the main menu, then press [ENTER]. A naming screen will be displayed:

```
ENTER: abcdefghi#
      Jklmnopqr#
Name:  #stuvwxyz#
      . , - +/ ():=#
```

Move the cursor to the wanted character, press [ENTER] key to confirm, this character will pop out under the *Name*, digits [0] to [9] and [-] could be entered by pressing the key on the panel directly. You could name the program a max eight-character word consisting of any combination of letters, numbers, and punctuation marks. If the length of name is less than eight-character, move the cursor to # and press [ENTER] or press [#] key on the panel directly to accept a file name.

If the name is already in use by a stored program, a screen saying "Name in use" will be display.

```
Enter: abcdefghi#
      Jklmnopqr#
Name:  #stuvwxyz#
in Use; . , - +/ ():=#
```

If this happens, press [ENTER], then enter a different name. After the name has been entered, a screen requesting selection of a temperature control method will be displayed:

```
Enter: a1
Control Method:
BLOCK sim-tube
```

Select a control method, then press [ENTER].

Attention please, under Block control method, the temperature shown on LCD represents only the temperature of the block. And under Sim-Tube control method, the temperature you see represents the temperature of the sample.

Note: Under Block control method, the best holding time is about 1 min. And under Sim-tube control method, the best holding time is about 30 seconds.

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```
Step # 1
94.0C Hold 0m 30s
```

If you want to get other running information, press [▼] key, and you could see: Press [▼] key once more to return to the former display.

```
RUN:    a1
Est remain: 1h01m20s
Step # 1 Tot Cyc 35
0h18m20s Segment # 2
```

If you want to review the gradient temperature when the instrument is running a gradient step (model 4196), press [▼] key once more and LCD will display as below:

```
RUN:    a1
Gradient 30C → 50C
Gradient View:
Col1 30C left →
```

If the program is finished normally, the LCD will show as below:
Press [CANCEL] return to the main menu

```
RUN:
Running finishes at
1h10m0s Press Key
```

Pause

At any time in running, the protocol could be paused by pressing [PAUSE].

```
<<Pause>> 94.0
Timer //// 0h19m40s
```

Press [PAUSE] again, the protocol will continue. You will find if the temperature is ramping while you press [PAUSE], the pause function will not work till the next temperature point is reached.

Stop

At any time in running, the protocol could be stopped by pressing [STOP]. And the screen will ask you:

```
Press<STOP>to abort,
<ENTER> to run again
```

Press [STOP] again to confirm stop protocol, or press [ENTER] to continue.

Auto Restart

If a power failure occurs when a protocol is running, the instrument will hold the protocol in memory. When power is restored, the protocol will begin running again at the point at which it was stopped, and a notice about the power interruption will be displayed. The LCD window will show like below:

```
History: Power Lost!
```

The message could be cleared by pressing [CANCEL] key, the protocol's runtime screen will immediately be displayed.

The screen will show:

```

RUN:
A1      a2
a3      a4
a5      a6
  
```

If you select "Yes", then the heated lid will work (holding at 105 °C) when running. You don't need to add oil in the tube.

```

RUN:      a1
Use heated lid?
YES      no
  
```

IMPORTANT: Attention please, BE SURE NOT TO ADD OIL IN THE WELL OF BLOCK AT ANY TIME.

If you select "No", then you must add oil in the tube, otherwise the sample in the tube will be condensed at the top of the tube.

```

RUN:      a1
Vessel:   0.2 TUBE
0.5 tube  plate
  
```

If you select Sim-Tube control mode in your protocol, the screen will show at this time:

After you select the tube or plate you use, the screen will show:

```

RUN:      a1
Volume(ul): 0
  
```

Input the volume of your sample, then the instrument will initiate running program. If you turn-on the heated-lid, the lid will firstly enter pre-heating stage. When the heated lid reaches 100 °C, the instrument begins running program. The heated lid continues heating until it reaches 105 °C.

Whenever the program finishes normally or abnormally, the heated lid will shut off automatically. In addition, you could select LID in main menu to set the temperature that the heated lid turn-off automatically.

```

RUN:      a1
Lid temperature : 20 °C
  
```

If you select turn-off the heated lid, we have to remind you once again that you must add oil in the tube, at this time the protocol will be launched without lid Heating display.

When running, you could see the normal running display from the LCD window:

```

RUN:      a1
Hold 40C  Blk=39.9C
0m30s    ;Step 2
  
```

If the step is a gradient step (only 4196 model), you could also see an additional gradient information on the normal running display:

```

RUN:      a1
Grad      30C      50C
Hold      40C      Blk=39.9C
0m30s    Step 2
  
```

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Step # 1
Option? no YES

* [Depending on the chosen model, on display will visualize 4196 or 4096]

The screen disappears within 10 seconds. If the self-test does not detect any problems, the Main Menu is displayed:

At this time, you could start to use the control panel and open the hot cover to put the PCR tubes into wells.


MAIN:	<4196>
RUN	enter
list	edit
files	lid

Note: Be sure not to touch the block or hot cover metal surface in case of being burnt

Note: If the hot cover turns on, the protocol will begin after the hot cover reaches 105 °C, which needs several minutes

Important: Please put 4 tubes (empty ones are also ok) at four corner of the block to balance the hot cover.

HOW TO CREATE A NEW PCR AMPLIFICATION FILE?

Main

Let's begin with the main menu. After self-test is finished, the screen shows; at any situation, press [STOP] key to return main menu.

Main menu has six functions – Run, Enter, List, Edit, File and Lid. Press [▶], [◀], [▼] or [▲] key to select function, press [ENTER] key to confirm

MAIN:	<4196>
RUN	enter
lis	edit
files	lid

Run

Under Main Menu, select Run and the screen will show all programs stored. Press [▼] or [▲] key to page-up or page-down, or press [▶] or [◀] key to move the cursor to look for the program you want to run, press [ENTER] key to confirm.

RUN:	
A1	a2
a3	a4
a5	a6

4. INSTALLATION / SETTING UP

Preliminary inspection

Unwrap the Thermal Cycler, take off the involving plastic and take off the polispan protection in which it comes fitted. Take off all the protective items and, without connecting the equipment to the net, make sure that it does not present any damage because of the shipment. In case the Thermal Cycler presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- Power cord
- Operation Manual
- Warranty Certificate

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.

Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT

Please put the Thermal Cycler on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the equipment near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.

Avoid inflammable or toxic substances in the working area.

The Thermal Cycler is supplied with a Schuko standard wire

Please insert the wire that feeds the AC electric current in the base of current 220 V, 50 Hz \pm 10% provided with earth wire and to the other end to the Thermal Cycler connector.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz \pm 10%.

Setting up

1. Turning the Thermal Cycler On

Move the power switch to "I" (the "On" position) which is located at the back of the machine, a self-test of the heat pump will usually begin running (see below).

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3. TECHNICAL SPECIFICATIONS


Code	50604096	50604196
Technology	Peltier	
Block		
Capacity	0.2 mLx96 wells + 0.5 mLx77 wells	
Temperature range	0° C ~ 99.9° C	
Heating rate	3° C/s	
Cooling rate	2° C/s	
Uniformity	< ± 0.2° C	
Accuracy	< ± 0.2° C	
Temp. control mode	Simulated tube or block	
Ramping rate adjust.	0.1° C ~ 3° C	
LCD screen	Large – size LCD	
Max.no.stor. programs	125	
Max.no. cycles	99	
Temp. increm-decrement	0.1 ~ 10° C	
Time increm-decrement	1 ~ 60 s	
Auto restart	Yes	
Hold at 4 °C	Always	
Auto pause	Yes	
Running time display	Yes	
Hot Cover		
Temperature	105° C	
Height	Adjustable	
Auto shut-off	Yes	
Gradient		
Accuracy	-	<± 0.3° C
Column uniformity	-	<± 0.3° C
Temperature range	-	30° C ~ 99.9° C
Temp. differential range	-	1° C ~ 30° C
Capacity	-	12 columns
Dimensions (LxWxH)	315x240x275 mm	
Weight	8.8 kg	
Power supply	110V±10%, 220V±25%, 50/60Hz	

* The manufacturer reserves the rights to change the specifications to improve the quality of products.

Control panel

- 2.1 [0] - [9] Digit key, for file parameter setting.
- 2.2 [STOP] Terminates a running protocol, and return to the main menu.
- 2.3 [PAUSE] Pause a running protocol. Press once more, the protocol will continue.
- 2.4 [CANCEL] Clear the present parameter.
- 2.5 [ENTER] Accepting the current setting, or operating the selection.
- 2.6 [▼] Page-down key.
- 2.7 [▲] Page-up key.
- 2.8 [◀] [▶] Cursor keys. Moving cursor from the present position to the right or the left.
- 2.9 [#] Name reception key. Press this key to finish the name of a file

```
LIST:      a1
Estimated run time:
99h99m99s
Save? yes NO
```

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```
Edit:      a1
Estimated Run Time:
1h 08m07s
Save?      YES no
```

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INDEX OF CONTENTS

1. USES OF THE INSTRUMENT	22
2. DESCRIPTION	23
3. TECHNICAL SPECIFICATIONS	25
4. INSTALLATION / SETTING UP.....	26
5. MAINTENANCE AND CLEANING	38
ANNEX I: CE CERTIFICATE	40

1. USES OF THE INSTRUMENT

The Thermal - Cyclor will meet your need for an easy-to-use, reliable, and a simple programming. It has been developed to provide all the fundamental features for the DNA amplification by PCR. With:

- Easy-to-learn programming interface for quick and painless programming
- Choice of simulated tube or block temperature control mode
- Space-saving design for easy setup and transportation
- Hot cover for oil-free operation
- Compact, light and noiseless

Software Functions

- File editing, saving and resuming
- Option of simulated tube or block mode
- Instantly displaying the information at any phase of program execution
- Pause of program execution
- Stop of program execution
- Auto-restart in case of power failure

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Thermal Cyler Nahita. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Thermal Cyler Nahita, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the equipment avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the Thermal Cyler to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the equipment.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the Thermal Cyler warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the Thermal Cyler warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.