AUXILAB S.L.

CE DECLARATION OF CONFORMITY

HORIZONTAL ELECTROFORESIS CELL of Auxilab,S.L for the Directive of machines (89/392/CEE modified) and the regulations adopted for its transposition

NAME OF THE MANUFACTURER / IMPORTER:
AUXILAB, S.L.

ADDRESS: **Polígono Morea Norte, 8
31191 Beriáin (Navarra)**

WE STATE THAT:

**HORIZONTAL ELECTROFORESIS CELL
53000005**

Are designed and manufactured according to:

- ◆ Directive 89/392/CEE, including the modifications and the national regulations that transpose them.
- ◆ Directive 73/23/CEE, modified over the electric security.
- ◆ Directive 89/336/CEE, modified over the electromagnetic compatibility.
- ◆ And that the following harmonized rules have applied (or part of them):
UNE 292-1/-2/-2/A1, UNE-EN 1050, UNE-EN 614-1, UNE-EN 1037, UNE-EN 1088, UNE-EN 547, UNE-EN 953, UNE-EN 294, UNE-EN 418, UNE-EN 894-1, UNE-EN 894-2, UNE-EN 954-1, UNE-EN 60204-1, UNE 61010-1/A2, UNE-EN61010-2-051.

BERIAIN 04th JUNE 2007

Signed by: **ALFONSO AINCIBURU SANZ**
DIRECTOR/MANAGER

Polígono Morea Norte, 8 31191 Beriain (Navarra) - Spain. Tel. 948 310 513 Fax 948 312 071
Internet: www.auxilab.es · Email: correo@auxilab.es

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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ÍNDICE DE IDIOMAS

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
7. MAINTENANCE AND CLEANING

To get the best results and a higher duration of the horizontal electrophoresis cell it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ NEVER autoclave or dry-heat sterilize the apparatus or components.
- ◆ Do not expose the apparatus or components to phenol, benzene, acetone, halogenated hydrocarbon solvents or undiluted laboratory alcohols.
- ◆ Avoid prolonged exposure of the apparatus or components to the UV light.
- ◆ In case of breakage or damage of the electrodes proceed as follows for their replacement:
 - Remove the screw and the rubber ring located under the banana connector of the tank.
 - Remove the broken or damaged electrode by pulling upward from the connector.
 - Insert the new electrode; do not forget the rubber ring.
 - Put the screw and the rubber ring again and tighten them to form a leak-free seal.

Cleaning

- ◆ Wash all components gently with water and non abrasive detergents, and rinse well in deionized water. *Note: be careful not to damage the electrode wire when cleaning the electrophoresis tank.*
- ◆ Wipe dry with a soft cloth or allow to air dry.
- ◆ Never use scourers or substances that can grate as they damage the electrophoresis cell and produce an early ageing of the equipment.

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

TROUBLE	COMENTS
Bromophenol blue dye turns yellow (pH change) during electrophoresis. Results are uninterpretable.	<ul style="list-style-type: none"> - Check the pH of the electrophoresis buffer. Be sure to use Tris base not Tris HCl. - Mix the buffer periodically during electrophoresis.
Samples leak underneath the gel upon loading	- The bottom of the wells were torn when the comb was removed. To minimize adhesion of the agarose to the comb, remove the comb as soon as the agarose is solidified.
Gel melts or become soft near sample wells	- This is due to the combination of pH drift and high temperatures. Reduce the electrophoretic voltage.
Pronounce "smiling" along one edge of the gel is observed (corresponding bands in different lanes migrate slower toward one edge of the gel)	- Gel was cast or electrophoresed out of level. Use the "bull's eye" level to verify that the apparatus is levelled prior to gel casting and electrophoresis.
S-Shaped lanes (anomalous migration front results in lanes that are not running at a uniform speed)	<ul style="list-style-type: none"> - Mix the buffer periodically during electrophoresis - Use a low conductivity / high buffering capacity buffer (eg. from TAE 1x or TBE 1x to TBE 0.5x) - Reduce the salt concentration of the sample
Bands are not straight lines or parallel to the top edge of the gel	<ul style="list-style-type: none"> - Verify that the wells are free of particles and bubbles before and after loading samples - Verify that the agarose is completely dissolved before casting the gel - Remove any particulate matter from the agarose before casting the gel - Be sure that bubbles are not trapped against the comb during gel casting
All bands appear as doublets (each band is represented twice within the same lane)	<ul style="list-style-type: none"> - Concentrate the sample and use a thin gel with a thin comb - Avoid movements during gel photography - Reduce voltage. Band doublets may result due to denaturalization from excess heat from electrophoresing gel at high voltage

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

5.3 Ethidium bromide

For visualization of double-stranded DNA after electrophoresis, it is necessary to stain the gel with a 0.5 g/mL solution of ethidium bromide in deionized water. The gel should be transferred to the solution and left to stain for a time depending on gel thickness (10-15 min for 3 mm-thick gels and up to 1 h for 10 mm-thick gels). To reduce fluorescent background, the gel can be destained in deionized water for 10-15 min (thin gels) or 30 min (thick gels).

Ethidium bromide can also be added directly to the agarose before casting the gel; so that electrophoresis is performed in the presence of ethidium bromide. The concentration of ethidium bromide in the gel must also be 0.5 g/mL.

It is recommended to prepare a concentrate stock of ethidium bromide that must be stored in the dark.

CHART 7. Concentrated stock of ethidium bromide, 20000x

Component	Quantity	Concentration
Ethidium bromide	100 mg	10 mg/mL
Deionized water	0.37 g	---

NOTE: ethidium bromide is a powerful mutagenic agent. Always wear gloves when handling solutions, gels or equipment that has been in contact with this substance.

6. TROUBLESHOOTING

TROUBLE

Bubbles do not appear no the electrodes when DC voltage is connected

COMENTS

- Verify that the DC power supply is operating properly
- Verify that the lid is properly placed on the apparatus
- Verify that power leads are in perfect conditions
- Verify that electrodes are in perfect conditions

"Flaming bands" (excessive fluorescence appearing as a trail above the band)

- Reduce the amount of the sample
- Reduce the amount of protein and/or glycerol of the sample

ter resolution for 0.1-3 Kb fragments. TBE buffer has a higher buffering capacity and a lower conductivity and therefore is better suited for electrophoresis at high voltage (>150 V). Moreover, this solution generates less heat at equivalent voltage and does not permit a significant pH drift.

Due to its low buffering capacity, TAE buffer requires recirculation systems or mixing periodically for full-length electrophoresis, particularly at higher voltages. TAE buffer provides better results for supercoiled DNA.

It is recommended to prepare a concentrated stock from which a 1x solution will be prepared when necessary.

CHART 4. TBE (Tris borate/EDTA) buffer, 10x

Component	Quantity	Concentration
Tris base	121.1 g	1 M
Boris acid, anhydrous	61.8 g	1 M
Na ₂ EDTA·2H ₂ O	7.4 g	20 mM
Deionized water	Up to 1 L	----

Note: This is a 10x concentrated solution. Dilute with deionized water prior to use. Final pH should be 8.3 at 25 °C.

CHART 5. TAE (Tris acetate/EDTA) buffer, 10x

Component	Quantity	Concentration
Tris base	48.4 g	400 mM
Na ₂ EDTA·2H ₂ O	7.4 g	20 mM
Sodium acetate, anhydrous	16.4 g	200 mM
Glacial acetic acid	17.0 mL	296 mM
Deionized water	Up to 1 L	----

Note: This is a 10x concentrated solution. Dilute with deionized water prior to use. Final pH should be 7.8 at 25 °C.

5.2 Loading buffer

Loading buffer must be added to samples before loading the gel. This solution must have enough glycerol to be denser than the electrophoresis buffer.

It is recommended to prepare a concentrated stock and immediately prior to loading add this stock to samples so that the final concentration at the mixture is 1x. Stock solution must be kept at 4 °C.

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her and the resolution lower. High voltages are suitable for selecting and analyzing purity of samples.

b. Low voltage: the heat produced will be lower and resolution higher but the time of electrophoresis will be longer.

- Check electrophoresis by following the migration of the bromophenol blue dye. Migration must be always towards the positive pole. **Note: Do not move the cell during electrophoresis. Small amounts of gases will be produced; make sure the equipment is located in a well aired area.**

- When electrophoresis finishes, switch the power supply off and unplug the power leads.

◆ After electrophoresis

- Remove the lid

- Take the tray with the gel out and carefully slide the gel out of the tray for staining with ethidium bromide and subsequent visualization with UV light. **Note: agarose gels tear easily if not properly supported.**

- Properly discard the electrophoresis buffer. Do not reuse the buffer.

- Rinse the tank with deionized water.

- Remove any residual agarose from the tray by rising with deionized water.

Security

◆ The horizontal electrophoresis cell must be used by previously qualified staff that knows how the equipment works thanks to the user manual.

◆ Put the electrophoresis cell on top of a horizontal, plane and stable table.

◆ Do not put the electrophoresis cell near any warm supply (burners, blowlamps...). Avoid vibrations.

◆ Horizontal electrophoresis cell is intended to use with a DC power supply. An inadequate manipulation may cause burns, hurts and damages to the electrical installation.

◆ Never connect the electrophoresis cell to the power supply before starting the experiment. The electrical current supplied by the power supply is transmitted to the cell through the leads and the lid assembly.

◆ Never use the electrophoresis cell without the lid.

◆ Always turn off the power supply before opening and removing the lid.

◆ Certain reagents indicated for use in this manual are of a hazardous nature (eg. ethidium bromide, acetic acid, boric acid). User should strictly observe the safety regulations indicated for such products.

◆ The equipments used in these procedures (power supplies, UV lamps, etc.) should be used following the manufacturer safety recommendations.

5. ELECTROPHORESIS BUFFERS AND OTHER SOLUTIONS

5.1 Electrophoresis buffers

For electrophoresis of gels of the same concentration and at the same voltage TAE buffer provides better resolution of fragments > 4 Kb in length while TBE buffer provides bet

Gel thickness (mm)	Comb thickness (mm)	Tooth width (mm)	Vol/well (μ L)	
10	2.0	40	720	
		20	360	
		15	270	
	1.5	7	94	
		5	67	
		3	37	
	5	2.0	40	320
			20	160
			15	120
1.5		7	42	
		5	30	
		3	12	
3		2.0	40	160
			20	80
			15	60
	1.5	7	21	
		5	15	
		3	6	

Note: volumes given are approximate. Low-percentage gels and low-melting-point agarose gels may have lower sample volumes.

- Add loading buffer to samples before loading the gel. Use an automatic pipette and disposable tips to prepare samples by combining appropriate quantities of loading buffer and sample. Mixtures can be prepared in microcentrifuge tubes or on a piece of Parafilm.

- Load the samples into each individual well using the pipetting device. Avoid bubble formation. *Note: caution must be taken at this point to avoid damaging the bottom and walls of the wells with the pipette tip.*

◆ Electrophoresis

- Cover the electrophoresis tank with the lid. *Note: The lid is especially designed to fit in a unique position thus assuring that the black lead is always connected to the black electrode and the red lead to the red electrode.*

- Plug the power leads to a power supply.

- Switch the power supply on and select the desired voltage. Small bubbles will rise from each electrode when the power supply is properly connected. *Note: Voltage will depend on thickness, length and percentage of the gel and type of electrophoresis buffer used:*

a. High voltage: the electrophoresis time will be shorter but the heat will be high

Tamaño gel (mm)	Grosor gel (mm)	Volumen agarosa (mL)	Volumen solución tampón (tanque) (mL)
60x60	5	18	560
	3	11	490
60x120, 120x60	5	36	550
	3	22	480
120x120	10	140	530
	5	72	500
	3	43	460

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- In a bottle or Erlenmeyer flask prepare the volume of liquid agarose required. Eg: to prepare a 120x120x10 mm and 1% agarose gel, add 1.4 g of agarose powder to 140 mL of electrophoresis buffer.

- Weight the flask.

- Heat the mixture in a hot plate or microwave until it boils and mix the flask occasionally until the agarose is completely dissolved in the buffer.

- Weight the flask again and adjust to the original weight by adding deionized water to compensate for evaporation.

- Apply a water spurt or introduce the flask into a water bath to cool the liquid agarose to 50-60 °C. **Note:** *casting gels with the agarose above 60 °C may cause the tray bottom to deform.*

- Pour the agarose at 60 °C on the gel tray fitted into the mould. Use a pipette tip to distribute the agarose evenly over the whole surface of the tray and to remove air bubbles especially those situated between and around the teeth of the comb.

- Allow the agarose to cool until completely solidified, usually 15-30 min.

- Remove the comb gently by pulling it vertically; take the tray with the gel out of the mould and put them at the central part of the electrophoresis tank. Make sure that the wells of the gel are near the cathode (negative pole, black) since the samples always migrate towards the anode (positive pole, red).

Note: *To store the gel, wet the gel surface with a small amount of electrophoresis buffer and seal the tray with the gel in a plastic bag. Keep at 4 °C, at these conditions the gel can be stored for 1 or 2 days.*

- Fill the tank with electrophoresis buffer until covering the surface of the gel with the volume indicated in Chart 2. **Note:** *the surface of the gel must be covered only with 1 or 2 mm of buffer. This avoids the gel to dry and assures an even voltage gradient through the gel. A higher volume of buffer is not necessary and will result in an increase of the current and heat.*

◆ Sample preparation and loading the gel

The amount of sample to be load is variable and depends on sample concentration and well size. Overloading the gel causes trailing and distortion of bands.

- Select the amount of sample to be loaded, bearing in mind the size of the wells (Chart 3) and the concentration of samples.

CHART 3. Maximum volume of sample depending on the type and thickness of the comb and the thickness of the gel.


◆ Set up

- Put the gel tray inside the mould for gel casting. Make sure that the black band for visualizing the wells stays between the grooves to hold the comb.

- Select the comb with the appropriate thickness and number of wells for the experiment. Usually the 1.0 mm-thick comb is used but when fragile gels are cast it is advisable to use the 1.5 mm-thick comb.

- Put the comb on the mould so that each extreme of the comb fits into one of the lateral grooves of the mould.

◆ Gel casting

- Select the agarose concentration of the gel according to the size of the DNA fragments to be analyzed (Chart 1).

CHART 1. Agarose concentration recommended for different sizes of DNA fragments.

Agarose concentration (w/v)	Size of DNA fragments (Kb)
0.5 %	1.0 to 30
0.7 %	0.8 to 12
1.0 %	0.5 to 10
1.2 %	0.4 to 3
1.5 %	0.2 to 3
2.0 %	0.01 to 1

- Select the volume of liquid agarose depending on the thickness of the gel to be cast (Chart 2).

CHART 2. Volume of agarose required according to gel thickness. It is also shown the volume of electrophoresis buffer needed to fill the tank in each case.

Gel size (mm)	Gel thickness (mm)	Agarose volume (mL)	Buffer volume (tank) (mL)
60x60	5	18	560
	3	11	490
60x120, 120x60	5	36	550
	3	22	480
120x120	10	140	530
	5	72	500
	3	43	460

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4. OPERATING MODE

4.1 Preliminary inspection

Unwrap the horizontal electrophoresis cell, take off all protections and make sure it does not present any damage because of the shipment. In case it presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- 1 Tank with electrodes
- 1 Lid with fixed leads
- 1 Mould for gel casting
- Set of gel trays: 2 trays for casting 60x60 mm gels, 1 tray for 120x60 mm gels, 1 tray for 60x120 mm gels and 1 tray for 120x120 mm gels.
- Set of combs: 1 comb 1 mm-thick for 11/25 wells, 1 comb 1.5 mm-thick for 8/18 wells, 1 comb 1.5 mm-thick for 6/13 wells and 1 comb 2 mm-thick for 1/2/3 wells.
- Warranty
- User's manual


We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.

4.2 Operating mode


Note: the equipment is exclusively designed for laboratory purposes; ambient temperature must be 0-40 °C and relative humidity must not exceed the 80%. Prior to use, put the electrophoresis cell on a plane, smooth and vibration-free surface.

- *Mould for gel casting*: it presents a special design to be used at both sides. Side A is for casting gels of 120x120 mm while side B is for casting gels of 60x60 mm, 120x60 mm and 60x120 mm. Each side presents lateral grooves to fit and hold the combs and in side A up to 2 combs can be held in case of having to analyze a large amount of samples at the same time.


Mould for gel casting

Combs

1. Electrophoresis tank
2. Replaceable electrodes
3. Antiskid legs
4. Device to facilitate lid removal
5. Gel tray
6. Lid
7. Slots for heat dispersion
8. Fixed leads

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INDEX OF CONTENTS

1. USES OF THE INSTRUMENT 17
 2. DESCRIPTION AND COMPONENTS 17
 3. TECHNICAL SPECIFICATIONS 19
 4. OPERATING MODE 19
 5. ELECTROPHORESIS BUFFERS AND OTHER SOLUTIONS..... 23
 6. TROUBLESHOOTING 25
 7. MAINTENANCE AND CLEANING..... 27
 ANNEX I: CE CERTIFICATE..... 28

1. USES OF THE INSTRUMENT

Nahita horizontal electrophoresis cell is specially designed for separating, identifying and preparing DNA and RNA molecules in agarose gels. The equipment is provided complete with all the necessary accessories for preparation and electrophoresis of 60x60, 60x120, 120x60 and 120x120 mm gels.

2. DESCRIPTION AND COMPONENTS


Components

Nahita electrophoresis cell is supplied together with all the necessary components for preparation and electrophoresis of agarose gels:

- *Electrophoresis tank*: made of a high quality material resistant to breakage, chemicals and high pressure and temperatures. It is a durable, transparent and well-sealed device. It presents 4 antiskid legs and 2 replaceable electrodes made of pure platinum and resistant to corrosion and high temperatures.

- *Lid*: it is also made of a high quality and resistant material. It presents slots for heat dispersion and two fixed power leads with banana connector to be plugged to a power supply.

- *Gel tray*: the equipment is provided with a set of gel trays of the adequate dimensions to cast gels of the different sizes previously indicated. All trays present an integrated rule and a dark band to facilitate well visualization and thus make sample loading easier.

Gel tray

- *Combs*: the equipment is supplied with a set of 1.0, 1.5 and 2.0 mm-thick combs that can be used on both sides depending on the gel size and the number of samples to be analyzed.

Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita horizontal electrophoresis cell. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita horizontal electrophoresis cell, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the electrophoresis cell avoiding sudden movements, knocks, free fall of heavy / sharp objects on it.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Accessories (including their loss) are not covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

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