**ESTEREOMICROSCOPIOS  
STEREOMICROSCOPES  
STÉRÉOMICROSCOPES**

**Zuzi**

**Serie / Series / Série 230**

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Thank you for choosing this equipment. We sincerely wish that you enjoy your stereomicroscope; we highly recommend looking after this equipment according to what is stated in this manual.

Zuzidevelops its products being remarkable in the ergonomics and security of its user. The incorrect use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages. Please read the point of Maintenance, where we expose the security notes.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

**TO GET THE BEST RESULTS AND HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.**

Please take in mind the following:

- ◆ This manual should be available for all the users of these equipments
- ◆ We should handle it carefully avoiding the sudden movements, knocks, or free fall of heavy objects, and its manipulation with sharp objects.
- ◆ Never dismantle the different pieces of the stereomicroscope to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty.
- ◆ To avoid fire o electric discharges, avoid humid environments.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your distributor or our ZUZI Technical Assistant Department by email [asistencia@auxilab.es](mailto:asistencia@auxilab.es)
- ◆ Please do not forget filling the Warranty Certificate and send it before 15 days after the date of purchase in order to obtain the two-year warranty.
- ◆ Neither lamps nor fuses are included in warranty.
- ◆ The manufacturer reserves the right to modify or improve the manual and/or the equipment.

**ATTENTION!!**

**IF AN EQUIPMENT IS NOT PROPERLY CLEAN AND DISINFECTED IT WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

**INDEX OF LANGUAGES**

Spanish .....	1-11
English .....	12-21
French .....	22-31

**ANNEX I: CE CERTIFICATE**

**AUXILAB S.L.**

**CE DECLARATION OF CONFORMITY**  
**STEREOMICROSCOPE ZUZI SERIES 230 of Auxilab,S.L for the Directive of Machines (89/392/CEE modified) and for the regulations adpoted for their transposition**

NAME OF THE MANUFACTURER/IMPORTER:  
**AUXILAB, S.L.**

ADDRESS: **Polígono Morea Norte, 8  
 31191 Beriáin (Navarra)**

WE STATE THAT:  
**STEREOMICROSCOPES ZUZI SERIES 230  
 CODES 50230030, 50230031, 50230037, 50230038**

Are designed and manufactured according to:

- ◆ Directive 89/392/CEE, including the modifications and the national regulations that transpose them.
- ◆ Directive 73/23/CEE modified over the electric security.
- ◆ Directive 89/336/CEE modified over the electromagnetic compatibility.

And that the following harmonized rules have been applied (or part of them):

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BERIAIN 8th July 2004

Signed by: ALFONSO AINCIBURU SANZ  
 DIRECTOR/MANAGER

Polígono Morea Norte, 8 31191 Beriain (Navarra) - Spain. Tel. 948 310 513 Fax 948 312 071  
 Internet: [www.auxilab.es](http://www.auxilab.es) · Email: [correo@auxilab.es](mailto:correo@auxilab.es)

**9. ACCESSORIES**

Different eyepieces and objectives can be attached depending on special needs. You can take photographs with model 231 if you have photography accessories.

<i>Objectives</i>		<i>Eyepieces</i>	
<i>Code</i>	<i>Description</i>	<i>Code</i>	<i>Description</i>
90230200	Objective 0.5x	90230146	Eyepiece WF 15x
90230202	Objective 2x	90230148	Eyepiece WF 20x
		90230175	Micrometric eyepiece
		90230191	Photographic eyepiece 2.5x
		90230193	Photographic eyepiece 5x

*Adapters*

<i>Code</i>	<i>Description</i>
90230011	Photographic adapter for series 230
90230419	Video adapter 1x

**INDEX OF CONTENTS**

1. USES OF THE EQUIPMENT ..... 13  
 2. DESCRIPTION..... 13  
 3. TECHNICAL SPECIFICATIONS ..... 14  
 4. INSTALLATION/SETTING UP..... 14  
 5. MAINTENANCE AND CLEANING ..... 15  
 6. TROUBLESHOOTING..... 17  
 7. CHOOSING THE OBJECTIVES AND EYEPIECES ..... 18  
 8. CAUSES OF A DEFECTIVE IMAGE ..... 19  
 9. ACCESSORIES ..... 20  
 ANNEX I: CE CERTIFICATE ..... 21

**1. USES OF THE EQUIPMENT**


These instruments are useful to observe enlarged images of small items. The image is not inverted but stereoscopic and maintains the sample's original colour when seen through the eyepiece. Zuzi equipments will help you in most difficult tasks related to control, research and documentation in the best work conditions.

As more remarkable applications, we stand out:

- Education: primary, secondary, practical and university education.
- Industry: electronics, watchmaking, textile, jewellery, paintings, food, equipments, casting, graphics, in every check, fixture and repair process, and quality control.
- Agriculture and cattle farming, minerals, biologics, arts restoration, philately, dentists, optics, research, graphology, mycology...

These instruments are very easy to use and they offer a high quality of the observed image. Besides, it is possible to obtain real documents by taking photographs.

**2. DESCRIPTION**

1. Incident light adjusting knob
2. Stage clips.
3. Stage.
4. Base.
5. Off / incident light / transmitted light/ incident and transmitted light
6. Incident light.
7. Zoom control.
8. Focusing adjusting knob.
9. Tension adjusting ring
10. Head.
11. Control of dioptric adjustment
12. Eyepieces.

### 3. TECHNICAL SPECIFICATIONS

**Base:** Iron made, very stable.

**Switch:** Off/ incident light/ reflected light/ reflected and transmitted.

**Light:** Incident and transmitted.

- Incident (6 fig. 1) 6V 15W halogen with intensity control
- Transmitted, halogen 6V 12W.

**Stage:** Frosted glass of 75 mm diameter for transmitted illumination and black/ white contrast of 75 mm diameter for incident illumination.

**Focusing adjusting knob:** Bilateral, by pinion and rack focusing system and security stop at the end of the haul.

**Zoom control:** Bilateral.

**Head:** 45° inclined and rotary 360° (models 237, 238). Adjustable interpupillary distance, dioptric adjustment in both ocular tubes.

- Binocular on models 230 and 237.
- Trinocular, with a tube for photographic attachment on models 231 and 238.

**Protective rubber cups:** To protect your eyes.

### 4. OPTICAL SPECIFICATIONS

**Eyepieces:** WF (wide field) 10X; field's diameter of 20 mm.

**Zoom:** It enables a continuous adjust of the magnified image without losing sharpness. Its zoom control changes continuously the magnifying factor from 0.75X to 3.4X (models 230, 231) and from 0.75X to 4.5X (models 237, 238). The working distance is exactly the same for all zoom range. It is operative at both sides, always being absolutely parafocal and centred.

**Optional:**

- WF 20X eyepiece - 5X photographic eyepiece
- WF 15X eyepiece - 0.5X objective
- 2.5X photographic eyepiece - 2.0X objective

#### Models 230, 231

ZOOM	OBJECTIVES	WF10X		WF15X		WF20X	
		magnifications	Field diam. (mm)	magnifications	Field diam. (mm)	magnifications	Field diameter (mm)
0.75 X ↓	1.0X	7.5x-3.4x	26 - 5.8	15x -68x	20 - 4.4	15x - 68x	16 - 3.5
	0.5X	3.75x-17x	52 - 11.6	5.6x -25.5x	40 - 8.82	7.5x - 34x	32 - 7
	1.5X	11.25x-51x	17.3 -3.19	16.8x -76.5x	13.3 - 2.94	22.5x 102x	10.6 - 2.35
3.4 X	2.0X	15x - 68x	13 - 2.9	22.5x -102x	10 - 2.2	30x - 136x	8 - 1.76

optional

- ◆ If the transmitted illumination lamp blows, replace it with a 6V 12W halogen lamp.
- ◆ To replace it, unscrew all the four screws located at the bottom of the base, loosen the lamp holder screws and change the lamp making sure you tighten it back.
- ◆ If necessary, replace the fuse doing as this: untwist the protective knob at the back of the base (*fuse*), replace the fuse and place it back in the same position twisting the protective knob back again.
- ◆ The stereomicroscope is supplied with a Schuko standard wire and it should be connected to a current wire provided with an earth wire; it should be handy to be disconnected in case of emergency.

#### Cleaning

- ◆ For cleaning the stainless steel steam parts use always a special cloth that does not give off fluff. Do not use any cloth that can produce damages to the equipment.
- ◆ For cleaning the different pieces of the equipment, we highly recommend to use neutral products that do not have a pH superior than 6 (soaped water without abrasives).

#### ATTENTION!!

IF AN EQUIPMENT IS NOT PROPERLY CLEAN AND DISINFECTED IT WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

### 8. CAUSES OF A DEFECTIVE IMAGE

In case of obtaining a deficient image, check the following:

- ◆ That the illumination is well done; that the luminous intensity is neither excessive nor too weak.
- ◆ That the optical system is clean, proceeding like this: turn and extract the eyepieces to see if there are little drops of dust or dirt and if so, clean them. You must never dismantle the head, but you can clean it carefully, blowing with a pear the accessible parts. If you own an optional objective and parasite images rotate at the same time you rotate the objective, clean it by using a dry brush trying to eliminate the dust. You should observe the front part with a magnifier or with an inverted ocular for a conscientious cleaning.

◆ In routine observations use the eyepieces with the lowest magnifications and the most powerful objectives. Eyepieces with large magnifications are reserved for special situations, having in mind they decrease definition and they do not get a better resolution.

## 7. CHOOSING OBJECTIVES AND EYEPIECES

*For a good functioning of the stereomicroscope it is necessary to follow the recommendations and the processes of use we give you.*

*Note:* All the process of use mentioned before will not have any value unless you keep a continued and careful maintenance.

### General recommendations

- ◆ Please follow the processes of use of this manual
- ◆ This manual should be available for all users of these equipments
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The maintenance necessary for these stereomicroscopes is the same that for other optical instruments. Always must be kept clean and covered to avoid dust.
- ◆ Never try to repair the equipment by yourself, since it could produce a defective use of the whole equipment or harm the users (burns, wounds) or cause damages to the electric installation. Also you can lose the warranty.
- ◆ In case of damage or doubt it can be solved getting into touch with the Zuzi Technical Service.
- ◆ For any manipulation or cleaning of the equipment please make sure it is disconnected from the net.
- ◆ Never dismantle the different pieces of the stereomicroscope to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty.
- ◆ Users must never take lens apart. If the outer lens surfaces were dirty, clean them with a smooth cloth (avoid fuzzy cloths) previously dampen in xilol or toluene. To eliminate dust on lens, blow them with a pear or clean them with a soft natural hair brush.
- ◆ In mechanical components use non corrosive lubricant, to avoid damages in paintings or terminations, paying special attention to avoid optical components contact. Always use aneutral detergent.
- ◆ Keep the stereomicroscope in a dry and fresh environment, out of humidity or even dust.
- ◆ If an incident illumination lamp blows, replace it with a halogen 6V 15W Zuzi lamp, avoiding contact with bare hands.

**IMPORTANT:** before replacing a bulb or a fuse, please make sure you have your equipment unplugged. Check you use a 6V 15W halogen lamp for incident illumination because any other kind of lamp will cause a failure in the equipment.

◆ To replace the incident illumination lamp, unscrew both the illumination system holder screws. Unscrew the two lamp holder screws and then, replace the lamp. Put the lamp holder back in its socket and this in the base. Make sure you have well adjusted every component.

Model	Zoom	Objectives	Working distance
230	0.75x	1.0x	90 mm
231	3.40x	0.5x	150 mm
237	0.70x	1.5x	49 mm
238	4.50x	2.0x	34 mm

### Models 237, 238

ZOOM	OBJECTIVES	WF10X		WF20X	
		magnifications	Field diam. (mm)	magnifications	Field diameter (mm)
0.7 X	1.0X	7x-4.5x	32.8-5.1	14x - 90x	16.4-2.6
↓	0.5X	3.5x-22.5x	6.5-10.2	7x - 4.5x	32.8-5.2
	1.5X	10.5x-67.5x	2.6-3.4	21x 135x	10.9-1.7
4.5 X	2.0X	14x - 90x	16.4-2.6	28x - 180x	8.2-1.3

opcional

## 5. INSTALLATION / SETTING UP

### Preliminary inspection

- ◆ Unwrap the stereomicroscope, take of the involving plastic and take of the polyspan protection in which it comes fitted, please make sure that it does not present any damage because of the shipment. If it does please communicate to your transport agent or supplier immediately, so that they can make the claims in the correct time limit.
- ◆ Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.
- ◆ Please check that all the accessories are enclosed with the equipment:
  - Stage of ground glass.
  - Black and white stage.
  - Eyepieces WF 10x.
  - Protective rubber cups.
  - 2 Fuses 0,5 A.
  - Plastic cover with "Ura Technic" and "Zuzi" printed on it.
  - Instruction manual.
  - Warranty certificate.

*We will not accept any return of equipments in return unless they come on their original wrapping.*

### Installation

- ◆ Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

**PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.**

◆ Please put the stereomicroscope on top of a horizontal, flat, stable table making a free space at least at 30 cm per side. Do not put the equipment near any warm supply (burners, blowlamp...), do not expose it directly to the sun, etc. During its functioning, do not use the equipment near fluids that can give off steams or produce explosive or inflammable mixtures.

- ◆ The stereo microscope is supplied with a Schuko standard wire.
- ◆ Please insert the wire that feeds the A/C in the base of current 220V  $\pm 10\%$  50-60Hz provided with earth wire.

*The manufacturer and the distributor will not assume any responsibility for the damages produced to the equipment during its installing or damages to persons suffered by the wrong use of the electric connection. The tension should be 220V  $\pm 10\%$  50-60Hz.*

- ◆ Always use 6V, 15W original Zuzi lamps for incident illumination. Please, do not use lamps of a higher power or you will probably over heat your stereomicroscope.
- ◆ For transmitted illumination, please use 6 V12W Zuzi fluorescent light.
- ◆ Cover the stereomicroscope with the included plastic bag when it is not being used to preserve it from dust laying on optical parts.
- ◆ Please keep the original wrapping to cover it when you are not using the equipment for a long time or for a safer transport in case it needs an overhaul.

#### Setting up

- ◆ Plug the stereomicroscope to a 220V  $\pm 10\%$ , 50-60 Hz socket ; it must be provided with an earth wire.
- ◆ Choose the most suitable light source (incident and/ or transmitted light) according to the sample you are going to observe.
- ◆ To adjust the incident illumination strength, use the adjusting knob in the base. These models have a ground glass stage for transmitted illumination and other for black/ white incident illumination.
- ◆ Watch through both eyepieces at the same time and adjust the interpupilar distance until both images become one stereoscopic image.
- ◆ Place the sample in the middle of the stage; slowly, turn the focusing control until you see a sharp image. By alternating the vision first with one eye closed, then with the other, you will appreciate an important focusing difference; to correct it, you will have to use the dioptric adjustment in each ocular tube.
- ◆ These models have an important advantage: magnifications can be continuously adapted. This way you can get a whole image of the sample's topography with few magnifications, and then you can observe with more detail the part you are more interested in. At the same time you can change the field's diameter from a larger diameter to a shorter one, without changing the objectives. The zoom range in these models goes from 0.75X to 3.4X (model 230,231) and 0.7X to 4.5X (model 237, 238).
- ◆ To replace the stage remove the one you were using and put the new one.
- ◆ Once you have finished working with the equipment we strongly recommend you to unplug it, as the electric components that its adapter has consume electricity. This can age the stereomicroscope too soon, cause damages due to the heating output, as well as an unnecessary waste of energy, etc.

#### Installation of the optional objectives

- ◆ If you need more accurate magnifications, you can attach optional objectives, as: 0.5x and 2x.
- ◆ Remove the zoom's protective knob and screw the objective required. Zoom usage is exactly the same. To know the magnifications you can obtain, see the chart on *Optical specifications*.

#### Taking photographs (models 231 and 238, optional)

- ◆ Place the photographic attachment code 90230011 (model 231) or 90230018 (model 238) over the head's tube and screw it on This adapter is composed by a photographic tube, a photographic eyepiece and a viewfinder, where the focusing eyepiece is located. The viewfinder ends in a T thread where you have to screw the corresponding T mount ring (not included) and then, over it, screw the photographic camera.
- ◆ Pull the observation/ photo rod located in the head to place it in *photo* position (outwards). This way the light is directed towards the photography tube, withdrawing from the left eyepiece.
- ◆ With the focusing ocular dioptric ring turn it until you see the cross sharply. The field to be photographed has to be on the photographic focusing eyepiece frame
- ◆ Take photographs by pressing the camera's release.

#### Security

- ◆ The stereomicroscope has to be handled by qualified staff for its correct use, so please keep this manual handy to help them learn how to use it properly.
- ◆ Please put the stereomicroscope on top of a horizontal, flat, stable table making a free space at least at 30 cm per side.
- ◆ Do not put the stereomicroscope near any warm supply (burners, blowlamp...), neither expose it directly to the sun; avoid vibrations and dust.
- ◆ During its functioning, do not use the equipment near fluids that can give off steams or produce explosive or inflammable mixtures.
- ◆ Never dismantle the different pieces of the stereomicroscope to repair it yourself, since it could produce a defective use of the whole equipment or harm the users (by burns, wounds) or cause damages to the electric installation.
- ◆ If you are not going to use the equipment during a long period of time please make sure that it is unplugged from the net.
- ◆ Made under European regulations.

## 6. CHOOSING OBJECTIVES AND EYEPIECES

- ◆ The image seen loses surface and sharpness as we increase magnifications. This increase must be done by changing to higher magnification objectives, and not by changing the eyepieces, as the latter only magnify the image that the objective provides. The more magnifications the eyepiece has, the bigger will be the decrease on sharpness, clarity and surface that the image can undergo.