

Viscosímetros serie 801 | **Viscosimeters serie 801**Marca|**Trademark: Nahita**

El viscosímetro rotacional analógico 801 de Nahita es un instrumento ideal para la determinación de la viscosidad media de diversos tipos de fluidos como aceites comestibles, lubricantes, pinturas, barnices y productos similares. Para ello, dispone de 4 velocidades de rotación y está equipado con 4 husillos que permiten al usuario seleccionar las condiciones más adecuadas según la viscosidad del líquido a analizar.

Gracias a su diseño, el viscosímetro 801 es un equipo estable, eficaz y cómodo de usar. La base presenta 3 patas roscadas para conseguir una perfecta nivelación del equipo mientras que la columna, maciza y con sistema de piñón y cremallera, y la pinza de sujeción del cabezal, con tornillo y guía, permiten colocar cómodamente el cabezal a la altura deseada y asegurando una buena estabilidad del mismo. Además, el dispositivo basculante para la sujeción del husillo facilita el trabajo y aumenta la eficacia del instrumento.

El equipo se suministra en un robusto maletín metálico para mantenerlo siempre bien recogido y en perfectas condiciones de uso.

Nahita analogical rotational viscosimeter 801 is ideal to determine viscosity of several types of fluids such as edible oils, lubricants, paints, varnishes and similar products. For this purpose, it presents 4 rotational speeds and is equipped with 4 rotors to allow user to select the more adequate conditions according to sample viscosity. Thanks to its design, the viscosimeter 801 is a stable, efficient and easy-to-use instrument. The base presents 3 levelling legs to get a perfect levelling of the equipment while the column, solid and with pinion & rack system, and the head holding platform, with screw and guide, allow placing the head at the desired height assuring a good stability. Moreover, the bascule rotor holding device facilitates operation and increases efficiency.

The equipment is supplied in a sturdy metallic case to keep it always in perfect conditions.



Referencia / Code	50801000
Rango medida / Measuring range	10-100000 mPas (cp)
Precisión / Accuracy	1%
Husillos / Rotors	nº 1, 2, 3, 4
Velocidad motor / Motor speed	6, 12, 30 ó 60 rpm
Tolerancia	± 5% (en fluidos Newtonianos)
Tolerance	± 5% (in Newtonian fluids)