

Bombas de membrana | Membrane vacuum pumps

Marca|Trademark:Nahita



Ref. / Code: 59001010



Ref. / Code: 59001025

Las bombas de membrana se pueden usar en un amplio rango de aplicaciones de vacío o para la eliminación de gases en laboratorio de química o física. Están totalmente exentas de aceite, no necesitan agua y, por lo tanto, no producen ningún tipo de aguas residuales. Son muy fáciles de usar, gracias a su bajo mantenimiento. Cuentan con manómetro metálico analógico.

The membrane pumps can be used in a wide range of vacuum applications or for the gas elimination in chemical or physics laboratories. Membrane pumps are totally oil-free, they need no water and therefore they do not produce any residual waters. They are very easy to use, thanks to their low maintenance requirements. Membrane pumps count with metallic analogy manometer.

Referencia / Code	59001010	59001025
Presión de vacío / Vacuum pressure	638 mm Hg (Torr)/0.085 MPa/850 mbar	600mm Hg (Torr)/0.08 MPa/800mbar
Caudal / Free air displacement	10 L/min - 0.6 m3/h	25 L/min - 1.5 m3/h
Consumo / Power	20 W	40 W
Peso / Weight	4 Kg	7 Kg
Dimensiones / Size	245x135x245 mm	280x140x220 mm
Alimentación / Power supply	220-240v 50hz	