

Alcohómetros | Alcoholometers

Marca|Trademark: Nahita

Los alcohómetros son un tipo de densímetros destinados a medir riquezas alcohólicas de líquidos, proporcionando los resultados en porcentaje de alcohol en volumen. Nahita le ofrece esta serie de alcohómetros con escala según Gay Lussac (% volumen de alcohol) además de un alcohómetro que presenta también la escala de Cartier (escala arbitraria en la que al agua le corresponden 0 °Cartier mientras que al etanol absoluto le corresponden 44 °Cartier).

Alcoholometers are a type of hydrometers designed to measure alcoholic richness of liquids giving results in percentage of alcohol in volume. Nahita offers this series of alcoholometers according to Gay Lussac (% alcohol volume) as well as an alcoholometer that also presents the Cartier scale (arbitrary scale in which 0 °Cartier correspond to pure water and 44 °Cartier correspond to absolute ethanol).



Referencia Code	Rango Range	Divisiones Divisions	Longitud Length	Diámetro Diameter
71304100	0-10 % vol.	1/10	330 mm	31 mm
71304110	10-20 % vol.	1/10	330 mm	31 mm
71304120	20-30 % vol.	1/10	380 mm	31 mm
71304130	30-40 % vol.	1/10	315 mm	30 mm
71304140	40-50 % vol.	1/10	315 mm	30 mm
71304150	50-60 % vol.	1/10	315 mm	28 mm
71304160	60-70 % vol.	1/10	380 mm	25 mm
71304170	70-80 % vol.	1/10	380 mm	25 mm
71304180	80-90 % vol.	1/10	380 mm	22 mm
71304190	90-100 % vol.	1/10	380 mm	23 mm
71304030	0-30 % vol.	1	210 mm	23 mm
71304050	0-50 % vol.	1	260 mm	16 mm
71304060	50-100 % vol.	1	280 mm	16 mm
71304990	0-100 % vol.	1	290 mm	15 mm
71306000	0-100 % vol.	1	290 mm	15 mm
	10-44 °Cartier	1/2		